





# PPH


## Disposable Prolapse Hemorrhoids Stapler and Accessories



Papiionaceous adjusting nut



Circular adjusting nut



Anorectal transfixation tool



Hook



Anorectal introducer



Anorectal dilator

Model	External Diameter	Internal Diameter	Number of Staples	Height of Staples
GTPH-32	32	22	32	4.1
GTPH-33	33	22	32	4.1
GTPH-34	34	24	34	4.1

### Characters

1. Reasonable number of staples conduce to reduce postoperative bleeding.
2. Imported titanium staples with higher intensity and tensility.
3. Integral shaping anvil to ensure the shape of staples.
4. Larger volume of cartridge to contain more tissue.
5. Adjust the thickness of cutting to the condition of tissue prolapsed.
6. Minimum gap between anvil and cartridge to avoid the damage to tissue caused by over pressure.
7. Two types of adjusting nut for operator's choice.
8. Selective cut accessories retain the normal anal cushion and mucous bridge, maintain maximum fine feeling and contraction of the anal function, can effectively prevent anas tomtotic stenosis.
9. Using a specially crafted anal mirrors to form a different window of an open-loop, only hemorrhoid mucosa is exposed, more targeted.

### Applicable Operation

Procedure for prolapse and hemorrhoids and selective removal of mucosa over the dentate line.




#### Manufacture Address:




# TST

## Disposable Prolapse Hemorrhoids Stapler and Accessories



Papilionaceous adjusting nut



Circular adjusting nut



One opening anoscope




Double openings anoscope



Three openings anoscope



Anorectal dilator



Hook

### Characters

1. Reasonable number of staples conduce to reduce postoperative bleeding.
2. Imported titanium staples with higher intensity and tensility.
3. Integral shaping anvil to ensure the shape of staples.
4. Larger volume of cartridge to contain more tissue.
5. Adjust the thickness of cutting to the condition of tissue prolapsed.
6. Minimum gap between anvil and cartridge to avoid the damage to tissue caused by over pressure.
7. Two types of adjusting nut for operator's choice.
8. Selective cut accessories retain the normal anal cushion and mucous bridge, maintain maximum fine feeling and contraction of the anal function, can effectively prevent anas tomtotic stenosis.
9. Using a specially crafted anal mirrors to form a different window of an open-loop, only hemorrhoid mucosa is exposed, more targeted.

### Applicable Operation

Procedure for prolapse and hemorrhoids and selective removal of mucosa over the dentate line.

Model	External Diameter	Internal Diameter	Number of Staples	Height of Staples
GTPH-32	32	22	32	4.1
GTPH-33	33	22	32	4.1
GTPH-34	34	24	34	4.1



#### Manufacture Address: